


High-Throughput Sequencing: A Brief Introduction

**James Pease
14 July 2014**

High throughput sequencing ("Next Generation Sequencing")



Many sequencing technologies

illumina®

454
SEQUENCING

AB Applied Biosystems

ion torrent
by life technologies™

PACIFIC BIOSCIENCES®



<http://www.youtube.com/watch?v=I99aKKHcxC4>



<http://youtu.be/v8p4ph2MAvI?t=52s>

Millions (billions) of short reads

[illegible]

Paired end sequencing

sequence each fragment
from both ends



Normally, data will come as two files
(one for each end, left/right)

DNA or RNA?

DNA

the whole genome

non-coding regions
(introns, etc.)

costs more

RNA

only transcripts

expression data

alternatively spliced
isoforms


costs less

RNA-Seq is actually cDNA-Seq

Reverse Transcriptase linked to Polymerase Chain Reaction (RT-PCR)


1

mRNA + primer 1 (BamHI)



Short Read Assembly

DNA Short Reads



Contiguous sequence “contig”

contigs



scaffold



chromosome




Mapping reads

**Alignment of DNA (or RNA) short reads
to chromosomes (or transcripts)
from a reference genome**




Microarray

Requires a pre-designed panel of known transcript sequences



	RNA-Seq	microarray
Annotate genes	yes	no
All transcripts	yes*	no
Differential Expression	yes	yes
Genotyping	yes	no*

Read Length and Error Rates inversely correlated



DNA-Seq Options

Whole-genome

RAD-Seq:

Restriction site Associated DNA Markers

Stranded/unstranded

RNA-Seq Options

Whole-transcriptome

Small RNAs

Stranded/unstranded

Normalized/non-normalized

Choosing a Sequencing Technology

- Read Length**
- Error Rates**
- Software Options**
- Cost (time and money)**

Researchers often spend months or years painstakingly designing a study and collecting data, and then throw their data into a computer and try to analyse it in minutes. The data deserve better. The quick and reckless approach to data analysis often fails to identify important aspects of the data. You should become friends with them! Conducting data analysis is like drinking a fine wine. It is important to swirl and sniff the wine, to unpack the complex bouquet and to appreciate the experience. Gulping the wine doesn't work.

Daniel Wright (2003) BJEP

Isolated DNA/RNA



fragmentation



sequencing

GATCCAGTACTGA

AGTAGTCATGCCG

TAGTAGTAGT

CGTACGTGTACC

CATGCGTACGTAG

TCATGCTAGTACGTA